

# ATLAS™ RIGID ARM

One control tightens or loosens all articulating joints to desired position.

90° angle for an extended reach

360° swivel

Choice arm lengths available:  
280mm, 400mm, or 550mm

Ergonomic knob to actuate quick connect mechanism.

**ACCESS**  
SURGICAL  
INNOVATIONS

**Engineered to hold an assortment of devices used in MIS surgical procedures. Designed to be more resilient to unintentional movement.**

# ATLAS™ Straight Column Arms

- AA-0350 Short Length  
280mm w/ Square QC
- AA-0150 Medium Length  
400mm w/ Square QC
- AA-0151 Medium Length  
400mm w/ Hex QC
- AA-0155 Medium Length  
400mm w/ Rack Clamp
- AA-0153 Medium Length  
400mm w/ Instrument Camp
- AA-0157 Medium Length  
400mm w/ Instrument Holder- Stainless Steel
- AA-0158 Medium Length  
400mm w/ Instrument Holder- Plastic
- AA-0450 Long Length  
550mm w/ Square QC



## X-Rigid Arm (extra hold)

- AA-0550 X-Rigid Arm  
400mm w/ Square QC

# ATLAS™ L-Column Arms

- AA-0300 Short Length  
280mm w/ Square QC
- AA-0100 Medium Length  
400mm w/ Square QC
- AA-0101 Medium Length  
400mm w/ Hex QC
- AA-0104 Medium Length  
400mm w/ Instrument holder
- AA-0107 Medium Length  
400mm w/ Instrument Clamp- Stainless Steel
- AA-0108 Medium Length  
400mm w/ Instrument Clamp- Plastic
- AA-0400 Long Length  
550mm w/ Square QC



## X-Rigid Arm (extra hold)

- AA-0500 X-Rigid Arm  
400mm w/ Square QC

## Quick Connect (QC) Mechanisms



AA-2041 Instrument Holder- Square QC Fitting  
Plastic Tip

AA-2051 Instrument Holder- Hex QC Fitting  
Plastic Tip



AA-2040 Instrument Holder- Square QC Fitting  
Stainless Steel Tip

AA-2050 Instrument Holder- Hex QC Fitting  
Stainless Steel Tip




AA-0120 AXS Instrument Clamp w/ Square Tip  
AA-0130 AXS instrument Clamp w/ Hex Tip



AA-0125 Rack Clamp w/ Square Tip  
AA-0135 Rack Clamp w/ Hex Tip

## Table Clamps



AA-0200 Table Clamp  
Aluminum




AA-0202 Table Clamp  
Rotatable  
Aluminum



AA-0204 Table Clamp  
Rotatable  
Stainless Steel

The quick connect feature safely loosens/locks in various instrument attachments to firmly hold in positions.




AA-0125



AA-2050



AA-2051



AA-0120



141 California Avenue Suite #102  
Half Moon Bay, California 94019  
[www.axcesssurgical.com](http://www.axcesssurgical.com)

phone: 650.563.9106  
fax: 650.563.9143

email: [orders@axcesssurgical.com](mailto:orders@axcesssurgical.com)  
questions: [info@axcesssurgical.com](mailto:info@axcesssurgical.com)  
customer service: 650.563.9231

